

compel[®]

Pivot Benching



compel[®]

Simple sleek solution.

Achieve a cohesive aesthetic from private offices to open floor plan benching with the versatile Pivit Collection. Available in the same four colors as Pivit Desks, Pivit Benching allows you to match or contrast your private offices with a highly functional and flexible benching system. Each 30" deep work surface is available in 48", 60" and 72" and supported by open frame end legs and shared center legs. Add mobile peds in Pivit Silver, frosted acrylic center privacy screens, and a trough power solution in both National and Chicago code. Increase functionality with our monitor arms and clamp on power.






Simple.

Pivit Benching is designed to deliver simplicity to benching applications. Simple to spec, simple to purchase, simple to install. Select your desired size, configuration, finish and accessories.

Worksurface Sizes



Finishes



American Cherry

Cafe

Grey Ash

Basic White

Accessories



Privacy Screen



Monitor Arm



Clamp On Power



Sleek.

Clean lines, smart accessories, and stylish finishes create a sleek and versatile aesthetic for Pivit Benching,



Beveled open frame leg in signature Pivit Silver



Mobile Pedestal in Pivit Silver with optional cushion



Low storage credenza in basic white with optional cushion



Solution.

Pivit Benching delivers a complete solution. The trough houses data and power for both National and Chicago code, and mounts for privacy screen. The privacy screen features a power and accessory pass-through for cables, monitor arms and additional power solutions. Pair with Pivit chassis and open frame desking for a cohesive solution throughout your entire office.



Panel End Desk In Grey Ash



Open Frame Desk In American Cherry



The Best of All Worlds

With a small business belief in customer service, a big business approach to global sourcing, and an aesthetic geared towards a design-minded market, Compel allows you to deliver a look and experience you're proud of within the price and timeline you need.

compel

compel[®]